









# **PiNK Plus Prestained Protein Ladder**

Catalog Number	Unit Size	Reactions
PM005-0500	500 ul	

- Broad range: 10-175 kDa
- Ready-to-use: supplied in a loading buffer for direct loading on gels
- Easy to identify: includes the  $\sim$ 10,  $\sim$ 40 and  $\sim$ 90 kDa reference bands coupled with an blue dye Sharp bands

Storage: Stable for up to 2 weeks at 25°C.

Stable for up to 3 months at 4°C.

For long term storage, store at -20°C.

## Description

The PiNK Plus Prestained Protein Ladder contains 11 proteins that resolve into sharp, tight bands in the range of 10-175 kDa. The PiNK Plus Prestained Protein Ladder allows you to monitor molecular weight separation during electrophoresis, estimate molecular weights of proteins of interest, and evaluate western transfer efficiency.

#### **Contents**

Approximately  $0.2\sim0.4$  mg/ml of each protein in the buffer (20 mM Tris-phosphate, pH 7.5 at 25°C), 2 % SDS, 1 mM Dithiothreitol, 4.8 M Urea, and 12 % (v/v) Glycerol).

### **Applications**

- Monitoring of protein migration during SDS-polyacrylamide gel electrophoresis.
- Monitoring of protein transfer onto membranes during Western blotting.
- Sizing of proteins on SDS-polyacrylamide gels and Western blots.

#### **Quality Control**

Tested in SDS-polyacrylamide gel electrophoresis and Western blotting














# — PROTOCOL —

## **Recommendations for Loading**

- 1. Thaw the ladder either at room temperature or at 37-40°C for a few minutes to dissolve precipitated solids. Do not boil.
- 2. Mix thoroughly to ensure the solution is homogeneous.
- 3. Load the following volumes of the ladder on SDS-polyacrylamide gel:
  - 5 ul per well for mini-gels, 2.5 ul per well for blots
  - 10 ul per well for large gels, 5 ul per well for blots



% of migration	Tris Glycine Gel					4-12% Bis Tris Gel		3-8% Tris Acetate
0 % —	8 %	10 %	12 %	15 %	4-20 %	MOPS	MES	TA
10 % — 20 % — 30 % — 50 % — 60 % — 80 % — 90 % — 100 %	175 130 95 70 62 51 42	175 130 95 70 62 51 42	175 130 95 70 62 51 42 29 22	175 95 130 95 70 62 1 42 29 22 14 10.5	175 130 95 70 62 51 42 29 22	165 115 80 61 52 42 32 24 17	155 120 82 62 55 45 34 26 16 12	160 1115 85 65 55 45 40 27 18



