

# **Blu11** Prestained Protein Ladder

03 JAN 2022

Catalog Number	Size	Range
PMB11-0500	500 μΙ	10 to 180 kDa

### **Storage Conditions**

Stable for up to 2 weeks at 25°C.

Stable for up to 3 months at 4°C.

Stable for up to 24 months at -20°C.

## Description

The Blu11 Prestained Protein Ladder is a combination of 11 pre-stained proteins with molecular weights from 10 to 180 kDa. The 11 recombinant proteins are covalently coupled with blue chromophore, and 2 reference bands at 25 kDa and 72 kDa serve as reference bands. The Blu11 Prestained Protein Ladder keeps track of the size and separation of proteins during SDS-polyacrylamide gel electrophoresis, approximating protein size and validating Western transfer efficiency on PVDF, nylon, or nitrocellulose membranes.

### Kit Content(s)

Blu11 prestained protein ladder	500μl x 1 vial	
---------------------------------	----------------	--

## Storage Buffer

Approximately  $0.1\sim0.5$  mg/ml of each protein in buffer (20 mM Trisphosphate pH 7.5 at 25°C), 2% SDS, 0.2 mM Dithiothreitol, 3.6 M Urea, and 15% (v/v) Glycerol.

### Required materials but not provided

- Vertical Electrophoresis system
- Power supplies
- Vortex or equivalent
- Microcentrifuge


### **Reaction Setup**

3  $\mu$ l or 5  $\mu$ l per loading for clear visualization during electrophoresis on 15-well or 10-well mini-gel, respectively. 2.5  $\mu$ l per well for general Western transferring.



### Guide for Molecular Weight Estimation (kDa)

Refer to the Blu11 Prestained Protein Ladder patterns in various electrophoresis conditions:



### **Important notes**

The molecular weight of each protein (kDa) was measured against an unstained protein ladder in every electrophoresis condition. Additional data should be considered for a more accurate adjustment.

